

Banco de Células do Rio de Janeiro

Data Sheet

PAGE 1/2

BCRJ Code:	0007	
Cell Line:	20.6.5	
Species:	Rattus norvegicus	
Vulgar Name:	Rat	
Tissue:	Blood	
Cell Type:	Hybridoma: B Lymphocyte	
Morphology:	Lymphoblast	
Growth Properties:	Suspension	
Products:	Immunoglobulin; monoclonal antibody; TcR murine: ab against; V-beta 2: ab against	
Biosafety:	1	
Addtional Info:	This hybridoma cell line secretes monoclonal antibody agaisnt murine TcR(V-beta 2 molecule). The origin of this cell line should be acknowledged in all relevant publications.	
Culture Medium:	RPMI-1640 medium modified to contain 2 mM L-glutamine, 1 mM sodium pyruvate, 4500 mg/L glucose and fetal bovine serum to a final concentration of 10%.	
Subculturing:	Cultures can be maintained by addition of fresh medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at 1 x 10e5 viable cells/mL. Maintain cultures at a cell concentration between 1 x 10e5 and 1 x 10e6 cells/mL. NOTE: Do not allow the cell concentration to exceed 1 x 10e6 cells/mL.	
Subculturing Medium Renewal:	Every 2 to 3 days	

0



Banco de Células do Rio de Janeiro

	Data Sheet	PAGE 2/2
Culture Conditions:	Atmosphere: air, 95%; carbon dioxide (CO2), 5% Temperature: 3	7°C
Cryopreservation:	95% FBS + 5% DMSO (Dimethyl sulfoxide)	
Thawing Frozen Cells:	SAFETY PRECAUTION: It is strongly recommended to always wear protective gloves, clothing, and a full-face mask when handling fr vials. Some vials may leak when submerged in liquid nitrogen, all nitrogen to slowly enter the vial. Upon thawing, the conversion of nitrogen back to its gas phase may cause the vial to explode or ej cap with significant force, creating flying debris. 1. Thaw the vial by gently agitating it in a 37°C water bath. To min contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 2 minutes). 2. Remove the vial from the water bath as soon as its contents and thawed and decontaminate it by dipping in or spraying with 70% ethanol. From this point, all operations must be performed under aseptic conditions. 3. For cells sensitive to DMSO, it is recommended to remove the cryoprotective agent immediately. Transfer the vial contents to a centrifuge tube containing 9.0 mL of complete culture medium a centrifuge at approximately 125 × g for 5 to 7 minutes. 4. Discard the supernatant and resuspend the cell pellet in the recommended complete medium (see specific batch information appropriate dilution ratio). 5. Incubate the culture under appropriate atmospheric and temp conditions (see "Culture Conditions" for this cell line). NOTE: It is important to avoid excessive alkalinity of the medium cell recovery. To minimize this risk, it is recommended to place the culture vessel containing the growth medium in the incubator fo 15 minutes before adding the vial contents. This allows the medi stabilize at its normal pH (7.0 to 7.6).	rozen owing of liquid ject its nimize ng re r strict a nd n for the perature during he r at least

Cellosaurus:

Depositors:

CVCL D659

f

0

Alberto Nobrega, Universidade Federal do Rio de Janeiro