

## Banco de Células do Rio de Janeiro

## **Data Sheet**

**PAGE 1/2** 

**BCRJ Code:** 0192

**Cell Line:** OFCOL II

**Species:** Mus musculus

**Vulgar Name:** Mouse; Balb/C

Tissue: **Bone Marrow** 

Cell Type: Osteoblast

Morphology: Fibroblast-Like

Disease: Normal

**Growth Properties:** Adherent

**Biosafety:** 1

Dulbecco's Modified Eagle's Medium (DMEM) modified to contain 2 **Culture Medium:** mM L-glutamine, 4500 mg/L glucose and 10% of fetal bovine serum.

Remove medium, and rinse with PBS without calcium and magnesium. Remove the solution and add an additional 1 to 2 mL of trypsin-EDTA solution. Allow the flask to sit at room temperature (or at 37°C) until the cells detach. Add fresh culture medium, aspirate and dispense into

new culture flasks. NOTE: For more information on enzymatic dissociation and subculturing of cell lines consult Chapter 12 in Culture of Animal Cells, a manual of Basic Technique by R. Ian Freshney, 6th

edition, published by Alan R. Liss, N.Y., 2010.

**Subculturing Medium** 

Renewal:

**Subculturing:** 

Every 2 to 3 days

**Subculturing Subcultivation** 

fb.com/bancodecelulas/

Ratio:

1:4 to 1:8 is recommended.





## Banco de Células do Rio de Janeiro

**Data Sheet** 

**PAGE 2/2** 

**Culture Conditions:** 

Atmosphere: air, 95%; carbon dioxide (CO2), 5% Temperature: 37°C

**Cryopreservation:** 

**Thawing Frozen Cells:** 

95% FBS + 5% DMSO (Dimethyl sulfoxide)

SAFETY PRECAUTION: It is strongly recommended to always wear protective gloves, clothing, and a full-face mask when handling frozen vials. Some vials may leak when submerged in liquid nitrogen, allowing nitrogen to slowly enter the vial. Upon thawing, the conversion of liquid nitrogen back to its gas phase may cause the vial to explode or eject its cap with significant force, creating flying debris.

- 1. Thaw the vial by gently agitating it in a 37°C water bath. To minimize contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 2 minutes).
- 2. Remove the vial from the water bath as soon as its contents are thawed and decontaminate it by dipping in or spraying with 70% ethanol. From this point, all operations must be performed under strict aseptic conditions.
- 3. For cells sensitive to DMSO, it is recommended to remove the cryoprotective agent immediately. Transfer the vial contents to a centrifuge tube containing 9.0 mL of complete culture medium and centrifuge at approximately 125 × g for 5 to 7 minutes.
- 4. Discard the supernatant and resuspend the cell pellet in the recommended complete medium (see specific batch information for the appropriate dilution ratio).
- 5. Incubate the culture under appropriate atmospheric and temperature conditions (see "Culture Conditions" for this cell line).

NOTE: It is important to avoid excessive alkalinity of the medium during cell recovery. To minimize this risk, it is recommended to place the culture vessel containing the growth medium in the incubator for at least 15 minutes before adding the vial contents. This allows the medium to stabilize at its normal pH (7.0 to 7.6).

**Depositors:** 

Alex Balduíno - Universidade Federal do Rio de Janeiro

Cellosaurus:

CVCL HA27

@bcrj\_apabcam



