

Banco de Células do Rio de Janeiro

Data Sheet

PAGE 1/2

BCRJ Code: 0207

Cell Line: QCPN

Species: Mus musculus

Vulgar Name: Mouse; Balb/C

Tissue: Spleen

Cell Type: Hybridoma: B Lymphocyte

Morphology: Lymphoblast

Growth Properties: Suspension

Female Sex:

Biosafety: 1

Iscove's Modified Dulbecco's Medium (IMDM) contains 2 mM L-**Culture Medium:** glutamine, 4500 mg/L glucose and fetal bovine serum to a final

concentration of 20%.

Cultures can be maintained by addition or replacement of fresh **Subculturing:** medium. Recommended cell density between 0,2 x 10e6 and 1,0 x

10e6 cells/mL.

Subculturing Medium

Renewal:

Every other day

Culture Conditions: Atmosphere: air, 95%; carbon dioxide (CO2), 5% Temperature: 37°C

Cryopreservation: 95% FBS + 5% DMSO (Dimethyl sulfoxide)



Thawing Frozen Cells:

Depositors:

Banco de Células do Rio de Janeiro

Data Sheet PAGE 2/2

SAFETY PRECAUTION: It is strongly recommended to always wear protective gloves, clothing, and a full-face mask when handling frozen vials. Some vials may leak when submerged in liquid nitrogen, allowing nitrogen to slowly enter the vial. Upon thawing, the conversion of liquid nitrogen back to its gas phase may cause the vial to explode or eject its cap with significant force, creating flying debris.

- 1. Thaw the vial by gently agitating it in a 37°C water bath. To minimize contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 2 minutes).
- 2. Remove the vial from the water bath as soon as its contents are thawed and decontaminate it by dipping in or spraying with 70% ethanol. From this point, all operations must be performed under strict aseptic conditions.
- 3. For cells sensitive to DMSO, it is recommended to remove the cryoprotective agent immediately. Transfer the vial contents to a centrifuge tube containing 9.0 mL of complete culture medium and centrifuge at approximately 125 × g for 5 to 7 minutes.
- 4. Discard the supernatant and resuspend the cell pellet in the recommended complete medium (see specific batch information for the appropriate dilution ratio).
- 5. Incubate the culture under appropriate atmospheric and temperature conditions (see "Culture Conditions" for this cell line).

NOTE: It is important to avoid excessive alkalinity of the medium during cell recovery. To minimize this risk, it is recommended to place the culture vessel containing the growth medium in the incubator for at least 15 minutes before adding the vial contents. This allows the medium to stabilize at its normal pH (7.0 to 7.6).

Andréa Gonçalves Trentini-Universidade Federal de Santa Catarina

Cellosaurus: CVCL 0R15



